

Decide and GO!

Hemoperfusion Platform* Ready when needed

intuitive

touch screen interface

PuriFi

VETRES 300 mL Device
Emergency and
Critical Care

maintenance free pressure monitoring

in-line transducers

easy

pre-assembled line kit

automatic

rinsing procedure auto-levelling bubble catcher

* Trademark for Punki by Cyrosonbonts pending

versatile


temperature management option

mobile

by battery back-up option

Veterinary Use Only





TECHNICAL DATA

Class	IIb
Dia a al munan	TE 44

Blood pump [5 to 450] mL/min; resolution: 5 mL/min

 $Syringe \ pump \ Supported \ syringes: 30 \ mL \ | \ Flow \ rate: [0.20] \ mL \ / \ h$

Sensitivity 0.1 mL | Accuracy: +/- 0.2 mL or 2%

Arterial pressure (DOME type) Range: [-250, +120] mmHg | Sensitivity: 1 mmHg

Displayed resolution 5mmHg

Pre-adsorber pressure (DOME type) Range: [-100, +700] mmHg | Sensitivity: 1 mmHg

Displayed resolution 5 mmHg

Auto-Leveling bubble detector Ultrasonic sensor

Sensitivity 0.1 mL (at the flow rate of 100 mL / min)

User interface 5.7" touch screen display

Membrane keyboard Containing 3 keys for the primary operation functions

Rotary selector To set the parameters

BLOOD TUBING

VET HP tubing set standard priming volume 90 mL

VET HP tubing set small priming volume 60 mL

Kalos heating priming volume 15 mL

ELECTRICAL DATA

Power supply / Frequency 115V (+/-10%), 60 Hz, 110VA

DIMENSIONS AND WEIGHT

Equipment height	365 mm
Equipment width	375 mm
Equipment depth	300 mm
Equipment weight	26.4 lbs
Cart height	930 mm
Cart base	440x500 mm
Maximum pole heigth	1970 mm
Trolley weight	79 lbs

ENVIRONMENTAL WORKING CONDITIONS

TRANSPORT / STORAGE CONDITIONS

Temperature	Range 59° F - 95° F	Range - 4° F - 122° F
Humidity	35% - 75%	35% - 85%
Atmospheric Pressure	Range 375 - 795 mmHg	Range 375 - 795 mmHg

Distributor

CytoSorbents Corporation